

## 兔抗 CDX2 多克隆抗体

中文名称: 兔抗 CDX2 多克隆抗体

英文名称: Anti-CDX2 rabbit polyclonal antibody

别名: CDX3; CDX-3; CDX2/AS

抗原: CDX2

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CDX2
<b>Immunogen:</b>	Synthetic peptide of human CDX2
<b>Full name:</b>	caudal type homeobox 2

<b>Synonyms :</b>	CDX3; CDX-3; CDX2/AS
<b>SwissProt:</b>	Q99626
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	34 kDa
<b>WB Positive control:</b>	231, HT-29, PC3 and hela cell
<b>WB Recommended dilution:</b>	500-2000

