

兔抗 CDX4 多克隆抗体

中文名称：兔抗 CDX4 多克隆抗体

英文名称：Anti-CDX4 rabbit polyclonal antibody

别名：caudal type homeobox 4

相关类别：一抗

抗原：CDX4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anterior-posterior patterning and hematopoiesis.
Applications:	ELISA, IHC
Name of antibody:	CDX4
Immunogen:	Synthetic peptide of human CDX4
Full name:	caudal type homeobox 4
SwissProt:	O14627
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer; Human colorectal cancer

IHC Recommend dilution:	20-100
-------------------------	--------



