

兔抗 CEACAM5 多克隆抗体

中文名称：兔抗 CEACAM5 多克隆抗体

英文名称：Anti-CEACAM5 rabbit polyclonal antibody

别名：CEA, CD66e

相关类别：一抗

储存：冷冻（-20℃）

抗原：CEACAM5

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Members of the CEACAM subfamily, including CEACAM 5, belong to the CEA gene family. For general information on the CEA gene family, see CEACAM1. Cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. Receptor for E.coli Dr adhesins. Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.
Applications:	ELISA, IHC
Name of antibody:	CEACAM5
Immunogen:	Synthetic peptide of human CEACAM5
Full name:	Carcinoembryonic antigen-related cell adhesion molecule

	e 5
Synonyms :	CEA, CD66e
SwissProt:	P06731
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human breast cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

