

兔抗 CEACAM21 多克隆抗体

中文名称：兔抗 CEACAM21 多克隆抗体

英文名称： Anti-CEACAM21 rabbit polyclonal antibody

别名： CEACAM3, R29124_1

抗原： CEACAM21

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	CEACAM21 is a gene which encodes carcinoembryonic antigen. It belongs to the immunoglobulin superfamily-CEA family. Three isoforms of the human protein are produced by alternative splicing. High levels of expression in hematopoietic system, brain and muscle.
Applications:	ELISA, IHC
Name of antibody:	CEACAM21
Immunogen:	Fusion protein of human CEACAM21
Full name:	carcinoembryonic antigen-related cell adhesion molecule 21
Synonyms :	CEACAM3, R29124_1

SwissProt:	Q3KPI0
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human tonsil
IHC Recommend dilution:	50-200

