

兔抗 EFTUD2 多克隆抗体

中文名称：兔抗 EFTUD2 多克隆抗体

英文名称： Anti-EFTUD2 rabbit polyclonal antibody

别名： elongation factor Tu GTP binding domain containing 2; MFDM; MFDGA; Snu114; Snrp116; SNRNP116; U5-116KD

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： EFTUD2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a GTPase which is a component of the spliceosome complex which processes precursor mRNAs to produce mature mRNAs. Mutations in this gene are associated with mandibulofacial dysostosis with microcephaly. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	EFTUD2
Immunogen:	Fusion protein of human EFTUD2

Full name:	elongation factor Tu GTP binding domain containin g 2
Synonyms:	MFDM; MFDGA; Snu114; Snrp116; SNRNP116; U5-1 16KD
SwissProt:	Q15029
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	109 kDa
WB Positive control:	Mouse brain tissue and Jurkat cell lysates
WB Recommended dilution:	500-2000



