

兔抗 EGFL7 多克隆抗体

中文名称：兔抗 EGFL7 多克隆抗体

英文名称： Anti-EGFL7 rabbit polyclonal antibody

别名： NEU1, ZNEU1, VE-STATIN, RP11-251M1.2

抗原： EGFL7

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a secreted endothelial cell protein that contains two epidermal growth factor-like domains. The encoded protein may play a role in regulating vasculogenesis. This protein may be involved in the growth and proliferation of tumor cells. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	EGFL7
Immunogen:	Fusion protein of human EGFL7
Full name:	EGF-like-domain, multiple 7
Synonyms :	NEU1, ZNEU1, VE-STATIN, RP11-251M1.2

SwissProt:	Q9UHF1
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

