

CDK2AP2 抗原（重组蛋白）

中文名称： CDK2AP2 抗原（重组蛋白）

英文名称： CDK2AP2 Antigen (Recombinant Protein)

别名： p14; DOC-1R

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	cyclin dependent kinase 2 associated protein 2
Synonyms:	p14; DOC-1R
Swissprot:	O75956
Gene Accession:	BC002850
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that interacts with cyclin-dependent kinase 2 associated protein 1. Pseudogenes associated with this gene are located on chromosomes 7 and 9. Alternatively spliced transcript variants have been observed for this gene.