

## CDKAL1 抗原（重组蛋白）

中文名称： CDKAL1 抗原（重组蛋白）

英文名称： CDKAL1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	CDK5 regulatory subunit associated protein 1 like 1
<b>Swissprot:</b>	Q5VV42
<b>Gene Accession:</b>	BC064145
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a member of the methyltransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.