

兔抗 EGR12 多克隆抗体

中文名称：兔抗 EGR1/2 多克隆抗体

英文名称：Anti-EGR1/2 rabbit polyclonal antibody

别名：TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268/AT591; CMT1D; CMT4E; KROX20

相关类别：一抗

储存：冷冻 (-20℃) 避光

抗原：EGR1/2

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Transcriptional regulator. Recognizes and binds to the DNA sequence 5'-CGCCCCGC-3'(EGR-site). Activates the transcription of target genes whose products are required for mitogenesis and differentiation.
Applications:	WB
Name of antibody:	EGR1/2
Immunogen:	Synthesized peptide derived from internal of human EGR1/2.
Full name:	early growth response 1/2
Synonyms :	TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268/A

	T591; CMT1D; CMT4E; KROX20
SwissProt:	P18146/P11161
WB Predicted band size:	58/50 kDa
WB Positive control:	Hela cells and Jurkat cells lysates
WB Recommended dilution:	500-3000

