

兔抗 EIF1AX 多克隆抗体

中文名称：兔抗 EIF1AX 多克隆抗体

英文名称： Anti-EIF1AX rabbit polyclonal antibody

别名： eukaryotic translation initiation factor 1A X-linked; EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： EIF1AX

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNA _i and eIF3) to the 5' end of capped RNA.
Applications:	ELISA, IHC
Name of antibody:	EIF1AX
Immunogen:	Fusion protein of human EIF1AX
Full name:	eukaryotic translation initiation factor 1A X-linked
Synonyms:	EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1
SwissProt:	P47813

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human cervical cancer
IHC Recommend dilution:	50-200



