

兔抗 EIF1AY 多克隆抗体

中文名称：兔抗 EIF1AY 多克隆抗体

英文名称： Anti-EIF1AY rabbit polyclonal antibody

别名： eukaryotic translation initiation factor 1A Y-linked; eIF-4C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： EIF1AY

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is located on the non-recombining region of the Y chromosome. It encodes a protein related to eukaryotic translation initiation factor 1A (EIF1A), which may function in stabilizing the binding of the initiator Met-tRNA to 40S ribosomal subunits. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	EIF1AY
Immunogen:	Fusion protein of human EIF1AY
Full name:	eukaryotic translation initiation factor 1A Y-linked
Synonyms:	eIF-4C

SwissProt:	O14602
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human liver cancer
IHC Recommend dilution:	50-200



