

DDA1 抗原（重组蛋白）

中文名称： DDA1 抗原（重组蛋白）

英文名称： DDA1 Antigen (Recombinant Protein)

别 名： PCIA1; C19orf58

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	DET1 and DDB1 associated 1
Synonyms:	PCIA1; C19orf58
Swissprot:	Q9BW61
Gene Accession:	BC000615
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins.