

FAM3A 抗原（重组蛋白）

中文名称： FAM3A 抗原（重组蛋白）

英文名称： FAM3A Antigen (Recombinant Protein)

别名： DLD, 2.19, XAP-7, DXS560S

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 34-230 amino acids of human FAM3A

技术规格

Full name:	family with sequence similarity 3, member A
Synonyms:	DLD, 2.19, XAP-7, DXS560S
Swissprot:	P98173
Gene Accession:	BC008912
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene belongs to the FAM3 family. It may act as a defensin against invading fungal microorganisms. Mostly expression in testis, pancreas, adrenal, placenta, brain, fetal brain, liver, kidney, skeletal muscle and heart.