

## LGR5 抗原（重组蛋白）

中文名称： LGR5 抗原（重组蛋白）

英文名称： LGR5 Antigen (Recombinant Protein)

别名： FEX, HG38, GPR49, GPR67, GRP49

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 22-320 amino acids of human LGR5

### 技术规格

<b>Full name:</b>	leucine-rich repeat containing G protein-coupled receptor 5
<b>Synonyms:</b>	FEX, HG38, GPR49, GPR67, GRP49
<b>Swissprot:</b>	O75473
<b>Gene Accession:</b>	BC096324
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes Orphan receptor. Stem cell marker of the intestinal epithelium and the hair follicle. Target gene of Wnt signaling.