

## 兔抗 CHAC2 多克隆抗体

中文名称：兔抗 CHAC2 多克隆抗体

英文名称：Anti-CHAC2 rabbit polyclonal antibody

别名：ChaC cation transport regulator homolog 2

抗原：CHAC2

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a gamma-glutamyl cyclotransferase that catalyzes the conversion of glutathione to 5-oxoproline and cysteinylglycine. It is thought that this gene is upregulated in response to endoplasmic reticulum stress and that the glutathione depletion enhances apoptosis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CHAC2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	ChaC cation transport regulator homolog 2
<b>SwissProt:</b>	Q8WUX2

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

