

兔抗 ATF5 多克隆抗体

中文名称：兔抗 ATF5 多克隆抗体

英文名称： Anti-ATF5 rabbit polyclonal antibody

别名： ATF5; HMFN0395

抗原： ATF5

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Transcriptional activator which binds the cAMP response element (CRE) (consensus: 5'-GTGACGT[AC][AG]-3'), a sequence present in many viral and cellular promoters and blocks the differentiation of neuroprogenitor cells into neurons. Its transcriptional activity is enhanced by CCND3 and slightly inhibited by CDK4.
Applications:	WB
Name of antibody:	ATF5
Immunogen:	Synthesized peptide derived from internal of human ATF5 .
Full name:	activating transcription factor 5

Synonyms :	ATFX; HMFN0395
SwissProt:	Q9Y2D1
WB Predicted band size:	31 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

