

MAP4K3 抗原（重组蛋白）

中文名称： MAP4K3 抗原（重组蛋白）

英文名称： MAP4K3 Antigen (Recombinant Protein)

别名： GLK; MEKKK3; MEKKK 3; MAPKKKK3; RAB8IPL1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 556-867 amino acids of human MAP4K3

技术规格

Full name:	Mitogen-activated protein kinase kinase kinase kinase 3
Synonyms:	GLK; MEKKK3; MEKKK 3; MAPKKKK3; RAB8IPL1
Swissprot:	Q8IVH8
Gene Accession:	BC071579
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the mitogen-activated protein kinase kinase kinase family. The encoded protein activates key effectors in cell signalling, among them c-Jun. Alternatively spliced transcripts encoding multiple isoforms have been observed for this gene. May play a role in the response to environmental stress. Appears to act upstream of the JUN N-terminal pathway. Ubiquitously expressed in all tissues examined, with high levels in heart, brain

n, placenta, skeletal muscle, kidney and pancreas and lower levels in lung and liver.