

MPZ 抗原（重组蛋白）

中文名称： MPZ 抗原（重组蛋白）

英文名称： MPZ Antigen (Recombinant Protein)

别名： PO; CHM; DSS; MPP; CMT1; CMT1B; CMT2I; CMT2J; CMT4E; CMTDI3; CMTDID; HMSNIB

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-248 amino acids of human MPZ

技术规格

Full name:	myelin protein zero
Synonyms:	PO; CHM; DSS; MPP; CMT1; CMT1B; CMT2I; CMT2J; CMT4E; CMTDI3; CMTDID; HMSNIB
Swissprot:	P25189
Gene Accession:	BC006491
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a major structural protein of peripheral myelin. Mutations in this gene result in the autosomal dominant form of Charcot-Marie-Tooth disease type 1 and other polyneuropathies. The myelin sheath is a multi-layered membrane, unique to the nervous system, that functions as an insulator to greatly

increase the velocity of axonal impulse conduction. The P0 glycoprotein, absent in the central nervous system,[1] is a major component of the myelin sheath in peripheral nerves. It comprises a large extracellular N-terminal domain, a single transmembrane (TM) region, and a smaller positively charged intracellular domain.