

## 兔抗 EMC1 多克隆抗体

中文名称：兔抗 EMC1 多克隆抗体

英文名称： Anti-EMC1 rabbit polyclonal antibody

别名： ER membrane protein complex subunit 1; CAVIPMR; KIAA0090

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： EMC1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a single-pass type I transmembrane protein, which is a subunit of the endoplasmic reticulum membrane protein complex (EMC). Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	EMC1
<b>Immunogen:</b>	Fusion protein of human EMC1
<b>Full name:</b>	ER membrane protein complex subunit 1
<b>Synonyms:</b>	CAVIPMR; KIAA0090

<b>SwissProt:</b>	Q8N766
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain and Human tonsil
<b>IHC Recommend dilution:</b>	40-200



