

兔抗 EMD 多克隆抗体

中文名称：兔抗 EMD 多克隆抗体

英文名称： Anti-EMD rabbit polyclonal antibody

别名： STA; EDMD; LEMD5

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： EMD

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Emerin is a serine-rich nuclear membrane protein and a member of the nuclear lamina-associated protein family. It mediates membrane anchorage to the cytoskeleton. D reifuss-Emery muscular dystrophy is an X-linked inherited degenerative myopathy resulting from mutation in the emerin gene.
Applications:	ELISA, WB
Name of antibody:	EMD
Immunogen:	Fusion protein of human EMD
Full name:	emerin
Synonyms :	STA; EDMD; LEMD5

SwissProt:	P50402
ELISA Recommended dilution:	2000-10000
WB Predicted band size:	29 kDa
WB Positive control:	293T and A549 cells
WB Recommended dilution:	1000-5000

