

## 兔抗 EME2 多克隆抗体

中文名称：兔抗 EME2 多克隆抗体

英文名称：Anti-EME2 rabbit polyclonal antibody

别名：essential meiotic structure-specific endonuclease subunit 2; SLX2B; gs125

相关类别：一抗

抗原：EME2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Interacts with MUS81 to form a DNA structure-specific endonuclease which cleaves substrates such as 3'-flap structures.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	EME2
<b>Immunogen:</b>	Synthetic peptide of human EME2
<b>Full name:</b>	essential meiotic structure-specific endonuclease subunit 2
<b>Synonyms:</b>	SLX2B; gs125
<b>SwissProt:</b>	A4GXA9
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human thyroid cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	231 cell lysate
<b>WB Recommended dilution:</b>	200-1000



