

兔抗 EMSY 多克隆抗体

中文名称：兔抗 EMSY 多克隆抗体

英文名称：Anti-EMSY rabbit polyclonal antibody

别名：EMSY, BRCA2 interacting transcriptional repressor; GL002; C11orf30

相关类别：一抗

抗原：EMSY

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Regulator which is able to repress transcription, possibly via its interaction with a multiprotein chromatin remodeling complex that modifies the chromatin. Its interaction with BRCA2 suggests that it may play a central role in the DNA repair function of BRCA2. As part of a histone H3-specific methyltransferase complex may mediate ligand-dependent transcriptional activation by nuclear hormone receptors.
Applications:	ELISA, IHC
Name of antibody:	EMSY
Immunogen:	Synthetic peptide of human EMSY

Full name:	EMSY, BRCA2 interacting transcriptional repressor
Synonyms:	GL002; C11orf30
SwissProt:	Q7Z589
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human esophagus cancer
IHC Recommend dilution:	30-150



