

MYLK2 抗原（重组蛋白）

中文名称： MYLK2 抗原（重组蛋白）

英文名称： MYLK2 Antigen (Recombinant Protein)

别名： KMLC; MLCK; MLCK2; skMLCK

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 297-596 amino acids of human MYLK2

技术规格

Full name:	myosin light chain kinase 2
Synonyms:	KMLC; MLCK; MLCK2; skMLCK
Swissprot:	Q9H1R3
Gene Accession:	BC069627
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a myosin light chain kinase, a calcium/calm odulin dependent enzyme, that is exclusively expressed in adult skeletal muscle.