

兔抗 EMX1 多克隆抗体

中文名称：兔抗 EMX1 多克隆抗体

英文名称：Anti-EMX1 rabbit polyclonal antibody

别名：empty spiracles homeobox 1

相关类别：一抗

储存：冷冻 (-20℃)

抗原：EMX1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Transcription factor, which in cooperation with EMX2, acts to generate the boundary between the roof and archipallium in the developing brain. May function in combinations with OTX1/2 to specify cell fates in the developing central nervous system.
Applications:	ELISA, IHC
Name of antibody:	EMX1
Immunogen:	Synthetic peptide of human EMX1
Full name:	empty spiracles homeobox 1
SwissProt:	Q04741
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human lung cancer and Human esophagus cancer
IHC Recommend dilution:	30-150



