

## NEUROG3 抗原（重组蛋白）

中文名称： NEUROG3 抗原（重组蛋白）

英文名称： NEUROG3 Antigen (Recombinant Protein)

别名： ngn3; Atoh5; NGN-3; Math4B; bHLHa7

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-214 amino acids of human NEUROG3

技术规格

<b>Full name:</b>	neurogenin 3
<b>Synonyms:</b>	ngn3; Atoh5; NGN-3; Math4B; bHLHa7
<b>Swissprot:</b>	Q9Y4Z2
<b>Gene Accession:</b>	BC117488
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a basic helix-loop-helix (bHLH) transcription factor involved in neurogenesis. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of congenital malabsorptive diarrhea 4 (DIAR4).