

## PLAU 抗原（重组蛋白）

中文名称： PLAU 抗原（重组蛋白）

英文名称： PLAU Antigen (Recombinant Protein)

别名： ATF; QPD; UPA; URK; u-PA; BDPLT5

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 179-417 amino acids of human PLAU

### 技术规格

<b>Full name:</b>	plasminogen activator, urokinase
<b>Synonyms:</b>	ATF; QPD; UPA; URK; u-PA; BDPLT5
<b>Swissprot:</b>	P00749
<b>Gene Accession:</b>	BC013575
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a serine protease involved in degradation of the extracellular matrix and possibly tumor cell migration and proliferation. A specific polymorphism in this gene may be associated with late-onset Alzheimer's disease and also with decreased affinity for fibrin-binding. This protein converts plasminogen to plasmin by specific cleavage of an Arg-Val bond in plasminogen. Plasmin in turn cleaves this protein at a Lys-Ile bond to form a two-chain derivative in w

hich a single disulfide bond connects the amino-terminal A-chain to the catalytically active, carboxy-terminal B-chain. This two-chain derivative is also called HMW-uPA (high molecular weight uPA). HMW-uPA can be further processed into LMW-uPA (low molecular weight uPA) by cleavage of chain A into a short chain A (A1) and an amino-terminal fragment. LMW-uPA is proteolytically active but does not bind to the uPA receptor. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.