

## PPP3R1 抗原（重组蛋白）

中文名称： PPP3R1 抗原（重组蛋白）

英文名称： PPP3R1 Antigen (Recombinant Protein)

别名： CNB; CNB1; CALNB1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	protein phosphatase 3 regulatory subunit B, alpha
<b>Synonyms:</b>	CNB; CNB1; CALNB1
<b>Swissprot:</b>	P63098
<b>Gene Accession:</b>	BC027913
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.