

## RETN 抗原（重组蛋白）

中文名称： RETN 抗原（重组蛋白）

英文名称： RETN Antigen (Recombinant Protein)

别名： ADSF, RSTN, XCP1, FIZZ3, RETN1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 19-108 amino acids of human RETN

技术规格

<b>Full name:</b>	resistin
<b>Synonyms:</b>	ADSF, RSTN, XCP1, FIZZ3, RETN1
<b>Swissprot:</b>	Q9HD89
<b>Gene Accession:</b>	BC069302
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene belongs to the family defined by the mouse resistin-like genes. The characteristic feature of this family is the C-terminal stretch of 10 cys residues with identical spacing. The mouse homolog of this protein is secreted by adipocytes, and may be the hormone potentially linking obesity to type II diabetes. Alternatively spliced transcript variants encoding the same protein have been found for this gene.