

SEMA3G 抗原（重组蛋白）

中文名称： SEMA3G 抗原（重组蛋白）

英文名称： SEMA3G Antigen (Recombinant Protein)

别名： sem2

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 533-782 amino acids of human SEMA3G

技术规格

Full name:	semaphorin 3G
Synonyms:	sem2
Swissprot:	Q9NS98
Gene Accession:	BC098137
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The transcription of this gene is activated by PPAR-gamma, and the resulting protein product plays a role in endothelial cell migration. Expression of this gene also inhibits tumor cell migration and invasion.