

SERPINE2 抗原（重组蛋白）

中文名称： SERPINE2 抗原（重组蛋白）

英文名称： SERPINE2 Antigen (Recombinant Protein)

别名： GDN; PI7; PN1; PNI; PI-7; PN-1; GDNPF

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 27-320 amino acids of human SERPINE2

技术规格

Full name:	serpin peptidase inhibitor, clade E (nexin, plasminogen activator inhibitor type 1), member 2
Synonyms:	GDN; PI7; PN1; PNI; PI-7; PN-1; GDNPF
Swissprot:	P07093
Gene Accession:	BC042628
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. Thrombin, urokinase, plasmin and trypsin are among the proteases that this family member can inhibit. This gene is a susceptibility gene for chronic obstructive pulmonary disease and for emphysema. Alternative splicing results in multiple transcript variants.