

## 兔抗 ERAP2 多克隆抗体

- 中文名称：兔抗 ERAP2 多克隆抗体
- 英文名称：Anti-ERAP2 rabbit polyclonal antibody
- 别名：endoplasmic reticulum aminopeptidase 2; LRAP; L-RAP
- 相关类别：一抗
- 抗原：ERAP2
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a zinc metalloaminopeptidase of the M1 protease family that resides in the endoplasmic reticulum and functions in N-terminal trimming antigenic epitopes for presentation by major histocompatibility complex (MHC) class I molecules. Certain mutations in this gene are associated with the inflammatory arthritis syndrome ankylosing spondylitis and pre-eclampsia. This gene is located adjacent to a closely related aminopeptidase gene on chromosome 5.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ERAP2

<b>Immunogen:</b>	Synthetic peptide of human ERAP2
<b>Full name:</b>	endoplasmic reticulum aminopeptidase 2
<b>Synonyms:</b>	LRAP; L-RAP
<b>SwissProt:</b>	Q6P179
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	10-50

