

ABI3BP 抗原（重组蛋白）

中文名称：ABI3BP 抗原（重组蛋白）

英文名称：ABI3BP Antigen (Recombinant Protein)

别名：TARSH; NESHBP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 731-930 amino acids of human ABI3BP

技术规格：

Full name:	ABI family, member 3 (NESH) binding protein
Synonyms:	TARSH; NESHBP
Swissprot:	Q7Z7G0
Gene Accession:	BC030221
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	A target of NESH-SH3/Abi3bp (TARSH) was originally identified as an SH3 domain-binding molecule of the NESH-SH3/Abi3 protein that is involved in Rac-dependent actin polymerization. In recent studies, TARSH gene expression was dramatically induced in mouse embryonic fibroblasts (MEFs) replicative senescence and suppressed in human lung carcinoma specimens and thyroid carcinomas. However, the molecular mechanism underlying the regulation of TARSH in tumorigenesis remains

ins unclear.