

## ACYP1 抗原（重组蛋白）

中文名称：ACYP1 抗原（重组蛋白）

英文名称：ACYP1 Antigen (Recombinant Protein)

别名：acylphosphatase 1; ACYPE

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 2-99 amino acids of human ACYP1

技术规格：

|                           |   |
|---------------------------|---|
| <b>Full name:</b>         | acylphosphatase 1   |
| <b>Synonyms:</b>          | ACYPE   |
| <b>Swissprot:</b>         | P07311  |
| <b>Gene Accession:</b>    | BC035568  |
| <b>Purity:</b>            | >85%, as determined by Coomassie blue stained SDS-PAGE  |
| <b>Expression system:</b> | Escherichia coli  |
| <b>Tags:</b>              | His tag C-Terminus, GST tag N-Terminus  |
| <b>Background:</b>        | This gene is a member of the acylphosphatase family. The encoded protein is a small cytosolic enzyme that catalyzes the hydrolysis of the carboxyl-phosphate bond of acylphosphates. Two isoenzymes have been isolated and described based on their tissue localization: erythrocyte (common) type acylphosphatase encoded by this gene, and muscle type acylphosphatase. Alternative splicing results in multiple transcript variants. |

**mlbio** 梅联生物  
Good elisakit producers