

ADAMDEC1 抗原（重组蛋白）

中文名称：ADAMDEC1 抗原（重组蛋白）

英文名称：ADAMDEC1 Antigen (Recombinant Protein)

别名：ADAM-like, decysin 1; M12.219

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 206-340 amino acids of human ADAMDEC1

技术规格：

Full name:	ADAM-like, decysin 1
Synonyms:	M12.219
Swissprot:	O15204
Gene Accession:	BC069582
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells.

mlbio 梅联生物
Good elisakit producers