

ADCK1 抗原（重组蛋白）

中文名称：ADCK1 抗原（重组蛋白）

英文名称：ADCK1 Antigen (Recombinant Protein)

别名：aarF domain containing kinase 1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 186-467 amino acids of human ADCK1

技术规格：

Full name:	aarF domain containing kinase 1
Swissprot:	Q86TW2
Gene Accession:	BC058906
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).