

AGFG1 抗原（重组蛋白）

中文名称：AGFG1 抗原（重组蛋白）

英文名称：AGFG1 Antigen (Recombinant Protein)

别名：ArfGAP with FG repeats 1; HRB; RAB; RIP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250amino acids of human AGFG1

技术规格：

Full name:	ArfGAP with FG repeats 1
Synonyms:	HRB; RAB; RIP
Swissprot:	P52594
Gene Accession:	BC096273
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene.

mlbio 梅联生物
Good elisakit producers