

AGTRAP 抗原（重组蛋白）

中文名称：AGTRAP 抗原（重组蛋白）

英文名称：AGTRAP Antigen (Recombinant Protein)

别名：angiotensin II receptor-associated protein; ATRAP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 77-159 amino acids of human AGTRAP

技术规格：

Full name:	angiotensin II receptor-associated protein
Synonyms:	ATRAP
Swissprot:	Q6RW13
Gene Accession:	BC017328
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates a ngiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.