

ALPI 抗原（重组蛋白）

中文名称：ALPI 抗原（重组蛋白）

英文名称：ALPI Antigen (Recombinant Protein)

别名：IAP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 203-503 amino acids of human ALPI

技术规格：

Full name:	alkaline phosphatase, intestinal
Synonyms:	IAP
Swissprot:	P09923
Gene Accession:	BC132678
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The intestinal alkaline phosphatase gene encodes a digestive brush-border enzyme. This enzyme is upregulated during small intestinal epithelial cell differentiation.