

ANGPTL7 抗原（重组蛋白）

中文名称：ANGPTL7 抗原（重组蛋白）

英文名称：ANGPTL7 Antigen (Recombinant Protein)

别名：angiopoietin-like 7; AngX; CDT6; dJ647M16.1; RP4-647M16.2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 122-343 amino acids of human ANGPTL7

技术规格：

Full name:	angiopoietin-like 7
Synonyms:	AngX; CDT6; dJ647M16.1; RP4-647M16.2
Swissprot:	O43827
Gene Accession:	BC001881
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Angiopoietin-related protein 7 is a protein that in humans is encoded by the ANGPTL7 gene, a member of the ANGPTL family. Highly and specifically expressed in the cornea where is confined to the stromal layer.

mlbio 梅联生物
Good elisakit producers