

ANKRD54 抗原（重组蛋白）

中文名称：ANKRD54 抗原（重组蛋白）

英文名称：ANKRD54 Antigen (Recombinant Protein)

别名：ankyrin repeat domain 54; LIAR

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 161 amino acids of human ANKRD54

技术规格：

Full name:	ankyrin repeat domain 54
Synonyms:	LIAR
Swissprot:	Q6NXT1
Gene Accession:	BC014641
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation. Interacts (via ankyrin repeat region) with LYN (via SH3-domain) in an activation-independent status of LYN. By similarity. Forms a multiprotein complex with LYN and HCLS1. By similarity. Interacts with TSN2, VAV1, DBNL AND LASP1.

mlbio 梅联生物
Good elisakit producers