

AOC1 抗原（重组蛋白）

中文名称：AOC1 抗原（重组蛋白）

英文名称：AOC1 Antigen (Recombinant Protein)

别名：amine oxidase, copper containing 1; ABP; DAO; KAO; ABP1; DAO1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 297 amino acids of human AOC1

技术规格：

Full name:	amine oxidase, copper containing 1
Synonyms:	ABP; DAO; KAO; ABP1; DAO1
Swissprot:	P19801
Gene Accession:	BC014093
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a metal-binding membrane glycoprotein that oxidatively deaminates putrescine, histamine, and related compounds. The encoded protein is inhibited by amiloride, a diuretic that acts by closing epithelial sodium ion channels. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

mlbio 梅联生物
Good elisakit producers