

## ARHGAP5 抗原（重组蛋白）

中文名称：ARHGAP5 抗原（重组蛋白）

英文名称：ARHGAP5 Antigen (Recombinant Protein)

别名：GFI2; p190-B; RhoGAP5; p190BRhoGAP

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 237 amino acids of human ARHGAP5

技术规格：

<b>Background:</b>	Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskeleton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isoforms have been found for this gene. GTPase-activating protein for Rho family members. May play a role in the reduction of the p21rasGTPase-activating potential of RASA1/p120GAP.
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Expression system:</b>	Escherichia coli
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Gene Accession:</b>	BC075799
<b>Swissprot:</b>	Q13017
<b>Synonyms:</b>	GFI2; p190-B; RhoGAP5; p190BRhoGAP
<b>Full name:</b>	Rho GTPase activating protein 5