

## ARHGAP15 抗原（重组蛋白）

中文名称：ARHGAP15 抗原（重组蛋白）

英文名称：ARHGAP15 Antigen (Recombinant Protein)

别名：BM046

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 280-470 amino acids of human ARHGAP15

技术规格：

<b>Full name:</b>	Rho GTPase activating protein 15
<b>Synonyms:</b>	BM046
<b>Swissprot:</b>	Q53QZ3
<b>Gene Accession:</b>	BC016701
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	RHO GTPases regulate diverse biologic processes, and their activity is regulated by RHO GTPase-activating proteins (GAPs), such as ARHGAP15. GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state. Has activity

toward RAC1. Overexpression results in an increase in actin stress fibers and cell contraction.