

ARHGEF5 抗原（重组蛋白）

中文名称：ARHGEF5 抗原（重组蛋白）

英文名称：ARHGEF5 Antigen (Recombinant Protein)

别名：P60; TIM; GEF5; TIM1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human ARHGEF5

技术规格：

Full name:	Rho guanine nucleotide exchange factor (GEF) 5
Synonyms:	P60; TIM; GEF5; TIM1
Swissprot:	Q12774
Gene Accession:	BC010046
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Rho GTPases play a fundamental role in numerous cellular processes initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form a complex with G proteins and stimulate Rho-dependent signals. This protein may be involved in the control of cytoskeletal organization.