

## ARL4A 抗原（重组蛋白）

中文名称：ARL4A 抗原（重组蛋白）

英文名称：ARL4A Antigen (Recombinant Protein)

别名：ARL4

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 2-197 amino acids of human ARL4A

技术规格：

<b>Full name:</b>	ADP-ribosylation factor-like 4A
<b>Synonyms:</b>	ARL4
<b>Swissprot:</b>	P40617
<b>Gene Accession:</b>	BC001111
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	ADP-ribosylation factor-like 4A is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL4A is similar to ARL4C and ARL4D and each has a nuclear localization signal and an unusually high guanine nucleotide exchange rate. ARL4A is located in both the nuclear and extranuclear cell compartments. Multiple transcript variants encoding the same protein have been found for this gene