

ARMCX3 抗原（重组蛋白）

中文名称：ARMCX3 抗原（重组蛋白）

英文名称：ARMCX3 Antigen (Recombinant Protein)

别名：ALEX3; dJ545K15.2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 111-379 amino acids of human ARM CX3

技术规格：

Full name:	armadillo repeat containing, X-linked 3
Synonyms:	ALEX3; dJ545K15.2
Swissprot:	Q9UH62
Gene Accession:	BC005194
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the ALEX family of protein s which may play a role in tumor suppression. The encoded protein contains a potential N-terminal transmembrane domain and a single Armadillo (arm) repeat. Other proteins containing the arm repeat are involved in development, maintenance of tissue integrity, and tumorigenesis. This gene is clos

ely localized with other family members on the X chromosome. Three transcript variants encoding the same protein have been identified for this gene.