

## ART3 抗原（重组蛋白）

中文名称：ART3 抗原（重组蛋白）

英文名称：ART3 Antigen (Recombinant Protein)

别名：ARTC3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 27-362 amino acids of human ART3

技术规格：

<b>Full name:</b>	ADP-ribosyltransferase 3
<b>Synonyms:</b>	ARTC3
<b>Swissprot:</b>	Q13508
<b>Gene Accession:</b>	BC008397
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes an arginine-specific ADP-ribosyltransferase. The encoded protein catalyzes a reversible reaction which modifies proteins by the addition or removal of ADP-ribose to an arginine residue to regulate the function of the modified protein. An ADP-ribosyl transferase pseudogene is located on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene.