

## ASB10 抗原（重组蛋白）

中文名称：ASB10 抗原（重组蛋白）

英文名称：ASB10 Antigen (Recombinant Protein)

别名：GLC1F

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 310-445 amino acids of human ASB10

技术规格：

<b>Full name:</b>	ankyrin repeat and SOCS box containing 10
<b>Synonyms:</b>	GLC1F
<b>Swissprot:</b>	Q8WXI3
<b>Gene Accession:</b>	BC126351
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. The SOCS box serves to couple suppressor of cytokine signaling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants have been described for this gene.