

## ASNS 抗原（重组蛋白）

中文名称：ASNS 抗原（重组蛋白）

英文名称：ASNS Antigen (Recombinant Protein)

别名：TS11

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human ASNS

技术规格：

<b>Full name:</b>	asparagine synthetase (glutamine-hydrolyzing)
<b>Synonyms:</b>	TS11
<b>Swissprot:</b>	P08243
<b>Gene Accession:</b>	BC014621
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is involved in the synthesis of asparagine. This gene complements a mutation in the temperature-sensitive hamster mutant ts11, which blocks progression through the G1 phase of the cell cycle at nonpermissive temperature. Alternatively spliced transcript variants have been described for this gene.