

## BCAP31 抗原(重组蛋白)

中文名称: BCAP31 抗原(重组蛋白)

英文名称: BCAP31 Antigen (Recombinant Protein)

别 名: CDM; BAP31; 6C6-AG; DXS1357E

储 存: 冷冻(-20℃)

相关类别: 抗原

概 述:

Fusion protein corresponding to a region derived from 124-246 amino acids of human BCAP31

## 技术规格:

Full name:	B-cell receptor-associated protein 31
Synonyms:	CDM; BAP31; 6C6-AG; DXS1357E
Swissprot:	P51572
Gene Accession:	BC014323
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the B-cell receptor associated prote in 31 superfamily. The encoded protein is a multi-pass transmembra ne protein of the endoplasmic reticulum that is involved in the anter ograde transport of membrane proteins from the endoplasmic reticul um to the Golgi and in caspase 8-mediated apoptosis. Microdeletion s in this gene are associated with contiguous ABCD1/DXS1375E delet ion syndrome (CADDS), a neonatal disorder. Alternative splicing of the



is	gene	results	s in	multip	ole :	transcript	variar	nts.	Two	related	pseudoge	n
es	have	been	ider	ntified	on	chromoso	ome 1	L6.				